

GENERATION OF TILED VIRUS AMPLICONS FOR RSV A SEQUENCING

1. TOPIC	2
2. STARTING MATERIAL	2
3. PROTOCOL.....	2
3.1. Generation of tiled amplicons	2
3.2. Purification of amplicons on beads	3
4. ANNEX 1	5

GENERATION OF TILED VIRUS AMPLICONS FOR RSV A SEQUENCING
1. TOPIC

Initial design by Etienne Simon-Loriere (etienne.simon-loriere@pasteur.fr) available at https://github.com/SimonLoriereLab/RSV_amplicons_panels

G5 Génomique Evolutive des Virus à ARN for Centre National de Référence – Virus des infections respiratoires (dont la grippe et le SARS-CoV-2), Institut Pasteur, France.

Adapted from: Quick, J. *et al.* Multiplex PCR method for MinION and Illumina sequencing of Zika and other virus genomes directly from clinical samples. *Nature Protocols* 12 (6), 1261-1276. <https://www.nature.com/nprot/journal/v12/n6/abs/nprot.2017.066.html>

Notes: The panel has been optimized using RSV A positive samples circulating in Europe in 2021-2022. Processing samples from previous years or from regions outside of Europe might require further adjustments.

2. STARTING MATERIAL

- RNA extracted with QiaAmp Viral RNA extraction kit (QIAGEN)
- Lunascript RT Supermix 5X kit (NEB)
- Q5 High-Fidelity PCR Kit (NEB) and dNTPs 10 mM
- Primers purchased at 200 µM from Sigma-Aldrich
- AMPure XP Bead-Based Reagent Protocol for PCR Purification (Beckman Coulter)

Equivalent kit and reagents can be used.

3. PROTOCOL
3.1. GENERATION OF TILED AMPLICONS

Reverse Transcription :

Composition	Reaction	Final Concentration
Lunascript Supermix 5X	2 µL	1X
RNA*	8 µL	-

*RNA was extracted from 100 µL of respiratory sample and eluted in 100 µL.

<u>LUNA RT</u>	
T°	Time
25°C	2 min
55°C	10 min
95°C	1 min
4°C	∞

GENERATION OF TILED VIRUS AMPLICONS FOR RSV A SEQUENCING

Amplicon generation by PCR in 2 independent reactions:

Primer pools.

Primers were purchased at 200 μ M from Sigma-Aldrich

Prepare two pre-mix (primers pools) with concentrations indicated in Annex 1

Composition	Reaction	Final concentration in 25 μ L
5X Q5 reaction buffer	5 μ L	1X
10mM dNTP	0,5 μ L	200 μ M
Primers pool 1 or primers pool 2 (see concentration in annex)	1,88 μ L of pre-mix pooled primers	See annex 1
cDNA	2,5 μ L	-
Q5 polymerase (2000 U/mL)	0,25 μ L	20 U/mL
H2O RNase DNase free	For final volume 25 μ L	-

Q5 PCR		
T°	Time	Nb Cycles
98°C	30 s	1
98°C	15 s	35
65°C	5 min	
4°C	∞	1

Check amplification on agarose gel:

Run 5 μ L of each reaction on a 1% agarose gel. Each should produce a visible ~400 bp band and the negative control (no cDNA) should be clean.

3.2. PURIFICATION OF AMPLICONS ON BEADS

Can use alternative purification techniques.

1. Vortex Agencourt AMPureXP beads bottle (placed at room temperature) to resuspend magnetic particles
2. For each sample, combine the PCR products from pool 1 and pool 2 PCRs (~50 μ L final volume)
3. Add 72 μ L (1,8X) of RNAClean XP beads to PCR; mix by pipetting 10 times
4. Incubate samples 10 minutes at room temperature.
5. Place the tubes on an appropriate magnetic rack to separate beads from the supernatant.
6. When the solution is clear (about 5-10 minutes), discard the supernatant without disturbing the beads.
7. Add 200 μ L of freshly prepared 70% ethanol to the tube while in the magnetic rack. Incubate at room temperature for 1 min, and then carefully remove and discard the supernatant
8. Repeat Step 7 once for a total of 2 washes.

GENERATION OF TILED VIRUS AMPLICONS FOR RSV A SEQUENCING

9. Completely remove the residual ethanol, and air-dry beads on the magnetic rack (5-10 minutes).
10. Add 40 μL of H₂O RNase Dnase free on the beads and remove the plate from the magnetic rack.
11. Mix well by pipetting up and down and put the tube in the magnetic rack until the solution is clear.
12. Transfer the 40 μL to a clean plate.
13. Quantify the DNA concentration using the Qubit High Sensitivity DNA kit (or equivalent) from 1 μL of each product. Expected range = 10-100 ng/ μL DNA.
14. Proceed with NGS library construction or other downstream application. Alternatively, the samples can be stored at -20°C .

The panel is compatible with Oxford Nanopore Technologies and Illumina sequencing. If using a combined paired read length shorter than the amplicons size (~400nt) with DNA fragmentation or tagmentation, adapt the bioinformatic pipeline for primer trimming. 150 nt long reads are the minimum recommended.

GENERATION OF TILED VIRUS AMPLICONS FOR RSV A SEQUENCING
4. ANNEX I

Name	Seq	Pool	Conc. Pre-Mix (μM)	Conc. Final (nM)
RSVA_1_LEFT	CTTGCGTAAACCAAAAAAATGGGG	1	0,4	30
RSVA_1_RIGHT	TGGTCATTGTTGAATCACTTAGTTTTTG	1	0,4	30
RSVA_3_LEFT	TAGAGACATCATAACACATAAATTTATATACTTG	1	0,2	15
RSVA_3_RIGHT	CAACTTGACTTTGCTAAGAGCCAT	1	0,2	15
RSVA_5_LEFT	AGGTATGTTATATGCTATGTCTAGATTAGGA	1	1	75
RSVA_5_RIGHT	TGGGTAATAAACCTTTATAGCGTTTCA	1	1	75
RSVA_7_LEFT	AAGTGATGTTACGGTGGGGG	1	0,2	15
RSVA_7_RIGHT	ACTCAAAGCTCTACATCATTATCTTTTGG	1	0,2	15
RSVA_9_LEFT	AAACCCAATAAATGAGACAGATGATACT	1	0,2	15
RSVA_9_RIGHT	TCCTGAGTCTTGCCATAGCTTC	1	0,2	15
RSVA_11_LEFT	GCTACCAAAGTGTACATCAAAACACA	1	0,1	8
RSVA_11_LEFT_bis	GCTACCAAAGTGTATATCAAAACACA	1	0,1	8
RSVA_11_RIGHT	TCTTAATGAGGGTCCCTTGGGT	1	0,2	15
RSVA_13_LEFT	ACACTCAACCCAACACATGACA	1	1,5	113
RSVA_13_LEFT_bis	ACACTCAACCCACACATGACA	1	0,5	38
RSVA_13_RIGHT	TTGTATGAATCTACTCATTGATGTAGAGG	1	1,5	113
RSVA_13_RIGHT_bis	TTGTGTGAATCTGCTCATTGATGTAGAGG	1	0,5	38
RSVA_15_LEFT	ACACGCTAGATAAAATCAACCAATGG	1	0,2	15
RSVA_15_RIGHT	CATTTGCCCCAGCGTTRTTYTT	1	0,2	15
RSVA_17_LEFT	CCACAAAGTCACACTAACAACCTGC	1	0,2	15
RSVA_17_RIGHT	TGGTGGTGGTTTTCTTTCCAGG	1	0,2	15
RSVA_19_LEFT	AAACCCTCCACTCAACCACCTC	1	0,2	15
RSVA_19_RIGHT	TCTGTACCATTACACTTATTTTCCTTGAT	1	0,2	15
RSVA_21_LEFT	GGCATTGCCGTATCCAAGGT	1	0,2	15
RSVA_21_RIGHT	TGCTGTCTAACTATTTGAACATTGCT	1	0,2	15
RSVA_23_LEFT	GTGACAATGCAGGATCAGTATC	1	0,2	15
RSVA_23_RIGHT	CAGAGGGGAACACTAATGGATC	1	0,2	15
RSVA_25_LEFT	TATACTGCAAGGCCAGAAGCAC	1	0,2	15
RSVA_25_RIGHT	CATGGGGTGGCCATTCAAATAA	1	0,2	15
RSVA_27_LEFT	ACTGAACTCAACAGCGATGACA	1	1	75
RSVA_27_RIGHT	CATTGAATGGTTGATCCGGTGG	1	1	75
RSVA_29_LEFT	CTCAGAATGTAACGCTTTAGGAAG	1	1	75
RSVA_29_RIGHT	GTAGTAATGACTGAGTTGTCTTCATC	1	1	75
RSVA_31_LEFT	GCAACATCCTCCATCATGGTTAATAC	1	1	75
RSVA_31_RIGHT	GTAGAACCCTCATTGTGGAATAG	1	1	75
RSVA_33_LEFT	TGCAGGTGACAATAACCTCAACA	1	0,4	30
RSVA_33_RIGHT	AAATTTTTAGGAGGTGATATAGCCTTATCA	1	0,4	30
RSVA_35_LEFT	AGGCAAAGAAAGAGAAGACTCAGTGT	1	0,2	15
RSVA_35_RIGHT	GGGGGTGCATGCCTATATGT	1	0,2	15
RSVA_37_LEFT	TCACTATTAGATCTAATATCTCTCAAAGGGA	1	2	150
RSVA_37_RIGHT	ACTGCATAATAGACTTTACCTCTATATTC	1	2	150
RSVA_39_LEFT	TCGAAGTTTCTATAGAAGAACTCCTGA	1	0,2	15

GENERATION OF TILED VIRUS AMPLICONS FOR RSV A SEQUENCING

RSVA_39_RIGHT	TCATAAACAACTCTTAGCCCGTGA	1	0,2	15
RSVA_41_LEFT	GAGATAAAAGAGAAATATTGAGTATGGAAAACC	1	1	75
RSVA_41_RIGHT	TGGCCTTCTCATATGTTAACCCA	1	1	75
RSVA_43_LEFT	TGGCCTTAGCTTAATGTCTGTAGT	1	0,2	15
RSVA_43_RIGHT	ACATATGATCAGTTATATATCCCTCTCCC	1	0,2	15
RSVA_45_LEFT	CATAGAGTAAAAGGATGTCATAGC	1	2	150
RSVA_45_RIGHT	CCTATAACAATAGTCACTCAGTGTC	1	2	150
RSVA_47_LEFT	CCTTAGTGCACAATAGCACATCAC	1	0,2	15
RSVA_47_RIGHT	CAGGCAATTCAGCATCACAGAC	1	0,2	15
RSVA_49_LEFT	GGTCCTGCAAATGTGTTCCCA	1	1	75
RSVA_49_RIGHT	AGCTGTTTAAACAATTCACCTTAGATGAGG	1	1	75
RSVA_2_LEFT	TGGCCTAATAGATGACAATTGTGAAA	2	0,4	30
RSVA_2_RIGHT	TGACCAGAAATGTAAATGTGGCCT	2	0,4	30
RSVA_4_LEFT	CAAACAATCCGAAATAACAACCTTTATGCA	2	0,2	15
RSVA_4_LEFT_bis	CAAACAATCCAAAATAACAGCTTTATGCA	2	0,2	15
RSVA_4_RIGHT	TCCACTCCATTTGCCTTAACATGA	2	0,2	15
RSVA_4_RIGHT_bis	TCCACTCCATTTGCTTTGACATGA	2	0,2	15
RSVA_6_LEFT	AGCAGCAGGAGATAGATCAGGT	2	0,2	15
RSVA_6_RIGHT	GGTTGTTCAATATATGGTAGAATCCTGC	2	0,2	15
RSVA_8_LEFT	AAATGGTGTGATTAACACTACAGTGATTAGA	2	0,15	11
RSVA_8_RIGHT	GGCGTAGGGTCTTCTTTGAAACT	2	0,15	11
RSVA_10_LEFT	AGGATTGATGAGAAATTAAGTGAATACTAGG	2	0,2	15
RSVA_10_RIGHT	TCAGCAGATAGATGTTTGGTTGGA	2	0,2	15
RSVA_12_LEFT	TCCATGCCAGCAGATCTACTCA	2	0,2	15
RSVA_12_RIGHT	TCAAGTGTGTTTCAGATCTTTATTTCTGA	2	0,2	15
RSVA_14_LEFT	CAACAAATTGGAAGCACACAGCT	2	0,2	15
RSVA_14_RIGHT	TGTTATCATGTGTATTAGTGTAAGTAAGGC	2	0,2	15
RSVA_16_LEFT	TCTCACCATGCAAGCCATYATC	2	0,2	15
RSVA_16_RIGHT	TGGGTGAGGTATGTTGGGGTT	2	0,2	15
RSVA_18_LEFT	ACCCAACAATGATTTTCACTTTGAAGT	2	0,2	15
RSVA_18_RIGHT	GTGTTGGATGAAGATRGAGRTTGTGAT	2	0,2	15
RSVA_20_LEFT	GCAAAGGCTATCTTAGTGCTCTAAG	2	0,2	15
RSVA_20_RIGHT	AGCCTTGTGTTGTGGATAGTAGAGC	2	0,2	15
RSVA_22_LEFT	GTGTAACACACCTGTAAGCAC	2	1	75
RSVA_22_RIGHT	CAGAGATTTACCTCACTTGGAATG	2	1	75
RSVA_24_LEFT	GGGTGGATACTGTGTCTGTAGG	2	0,2	15
RSVA_24_RIGHT	GGTGCTATTTTTATTTCAGTTACTAAATGCA	2	0,2	15
RSVA_26_LEFT	GGGGCAAATATGTCACGAAGGA	2	0,2	15
RSVA_26_RIGHT	GCTTCAATATATGATATGACAGTATTGTACAC	2	0,2	15
RSVA_28_LEFT	CAAGTAATTGTAGAGTCACTATGTATAACCA	2	0,2	15
RSVA_28_RIGHT	AGTTGGTATAATCATTTTTGAGATAAGGACC	2	0,2	15
RSVA_30_LEFT	CTTATGACATACAAGAGTATGACCTC	2	0,2	15
RSVA_30_RIGHT	CTCACTAGATCGATATTGTGTTAATATGC	2	0,2	15
RSVA_32_LEFT	TGTGGATTCAATAATGTTATCTTGACACA	2	0,2	15
RSVA_32_RIGHT	AGCATCCATGGCTTGCTTTCA	2	0,2	15
RSVA_34_LEFT	TCAGGACTACGTTTCTATCGAGAGT	2	0,2	15

GENERATION OF TILED VIRUS AMPLICONS FOR RSV A SEQUENCING

RSVA_34_RIGHT	TGTTTTCTGCTATCATTTTCTCTGCT	2	0,4	30
RSVA_36_LEFT	TGGATGAACTGCATGGTGTACA	2	0,4	30
RSVA_36_RIGHT	TGTGGCCTATTCCTGCATACTC	2	0,75	56
RSVA_38_LEFT	GGACCGTGGATAAACACTATACTTGA	2	0,25	19
RSVA_38_LEFT_bis	GGACCGTGGATAAACACCATACTTGA	2	1	75
RSVA_38_RIGHT	CTTGAAGTTTGTCTTTTAAATCATGGTTTG	2	0,2	15
RSVA_40_LEFT	GTGCACAACACTATAACCACCAC	2	0,2	15
RSVA_40_RIGHT	TGTCCATTGTATACATGATACTGGGT	2	0,2	15
RSVA_42_LEFT	AGCAAAATTGGATTGGGTGTATGC	2	0,2	15
RSVA_42_RIGHT	ACATCACCTGTGAATATGGGAGG	2	0,4	30
RSVA_44_LEFT	ACTAATTTAGCTGGACATTGGATTCTT	2	0,4	30
RSVA_44_RIGHT	TGTTAACAACCCAAGGGCAAAC	2	2	150
RSVA_46_LEFT	ACAACAAATTATATCATCCYACACCAG	2	2	150
RSVA_46_RIGHT	TATGATGCCAAGGAAGCATGC	2	0,4	30
RSVA_48_LEFT	GACCATTCTGCTACAGATGCA	2	0,4	30
RSVA_48_RIGHT	GGCATGATGAAATTTTTGGTTCTTGA	2	1	75
RSVA_50_LEFT	CGATAAAGAGTCTATTGATGCAAATATTAAGAG	2	1	75
RSVA_50_RIGHT	TGGATTTAACTTATTCTTCCTAGATCAAAA	2	0,4	30